



SOPRO 595

Specifications

Sopro 595

- CCD 1/3" high sensitivity
- Resolution: (752x582) PAL; (768x494) NTSC
- Definition: 470 lines
- Sensitivity: 2 lux
- Lighting: 8 LED
- Adjustments: automatic
- Non-inverted image
- Angle of view: 80°
- Cable length: 2.5m (option: 5m)
- Size of the handpiece: L: 198 x W: 27 x H: 18mm
- Size of the distal part: W: 16.5 x H: 12mm
- Weight: 50g

Dock 595M

- Storage of one or four images
- Power supply: 115V-60Hz & 230 V-50Hz
- Power consumption: 16VA
- One video output PAL or NTSC
- One S-video output PAL or NTSC
- Size of the dock: D: 190 x W: 100 x H: 40mm
- Weight of the dock: 400g

Dock 595MU

- Storage of one or four images
- Power supply: 24 VAC; 50-60Hz
- Power consumption: 16 VA
- One video output PAL or NTSC
- One S-video output PAL or NTSC
- Size of the dock: D: 135 x W: 80 x H: 35mm
- Weight of the dock: 250g

Dock 595

- Power supply: 115V-60Hz & 230 V-50Hz
- Power consumption: 12VA
- One video output PAL or NTSC
- One S-video output PAL or NTSC
- Size of the dock: D: 66 x W: 66 x H: 28 mm
- Weight of the dock: 70g

Dock 595USB

- One USB output
- Size of the dock: D: 75 x W: 29 x H: 17mm
- Weight of the dock: 86g

