


# SOPRO 717

## Specifications

### SOPRO 717

- High sensitivity 1/4" CCD
- Resolution: (752x582) PAL; (768x494) NTSC
- Definition: 470 lines
- Sensitivity: 2 lux
- Lighting: 8 LED
- Adjustment: 3 pre-set positions (Extra-oral, Intra-oral, Macro)
- 2 **SoproShade** positions
- Non-inverted image
- Freeze Frame with **SoproTouch**  or pedal (option)
- Angle of view: 70°
- Cable length: 2,5m

- Size of the handpiece:  
L.200 mm x l.28 mm x H.24mm
- Size of the distal part: l.13mm x H. 8mm
- Weight: 75g

### DOCKING STATION

- Storage of one or four images
- Power supply: 115V-60Hz & 230V-50Hz
- Power consumption: 9VA
- One PAL or NTSC video output
- One PAL or NTSC S-video output
- One digital USB2 output
- Size of the dock:  
L.145mm x l. 130mm x H. 35mm
- Weight of the dock: 245g

